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DE

Application and safety instructions for the MEISINGER

Implant Preparation Kit with Stop

Initial Bur System with Stop for safe Implant Site Preparation

Art.-No.: IPK02

The Implant Preparation Kit provides the easy and safe preparation of the implant site. The stop sleeves allow for a defined depth stop during the initial and pilot drilling. For the application without stop sleeves, the instruments are equipped with laser markings.

Indications for Use

The Bone Management® system Implant Preparation Kit with Stop is indicated for marking the implant site, performing the subsequent pilot drilling using depth-defining stop sleeves and the system independent preparation for surgical implant site preparation.

Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Determination of the drilling depth / selection of the stop sleeve

The required drilling depth and selection of the appropriate stop sleeve is determined before the start of the operation according to the selected implant length. With a required drilling depth of 8.0 mm for example, the stop sleeve CL038 must be selected. This has a length of 11.0 mm and allows a maximum drilling depth of 8.0 mm when mounted.

ATTENTION: In order to avoid injury to anatomical structures, the apical extra length of the pilot drill of 0.8 mm must be considered.

Preparation of the implant site

For the following steps, the selected stop sleeve must be placed on the respective drill used, as shown in the illustration. Make sure that the stop sleeve is securely seated.

To avoid overheating of the jawbone, ensure sufficient cooling with preferably physiological, sterile saline solution while using the drills. In addition, the drills are to be used intermittently, i.e. with up and down movements.

Pilot drilling of the alveolar ridge

In order to prevent the following instruments from slipping when, for example, the alveolar ridge tapers to a point, the pilot drill can be used beforehand.

Performing the pilot drilling

The pilot drilling is now carried out at the pre-drilled implantation site. This is used to create an implant site with a diameter of 2.0 mm. In addition to the predefined depth, which is specified by the attached stop sleeve, the alignment of the implant site is determined with the pilot drilling.

Implant site preparation and implant insertion

Subsequently, the extension to the planned implant diameter and the fine preparation of the implant site are carried out with the system's own instruments before the insertion of the selected implant can take place in compliance with the respective manufacturer's information.

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered. In addition, the contraindications for the implantation of an endosseous dental implant should be considered. The respective implant manufacturer provides information about these in the user instructions of the selected implant system. There are no specific contraindications for the use of the instruments contained in this kit.

IMPORTANT: Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent teeth / tooth roots (risk of damage, infection/dehiscesences)

CAUTION: The surgical and implantological procedure must always be carried out as per the current state of the art. For all drilling procedures, the following applies: In order to reduce the risk of overheating the bone and necrosis formation along with it, the respective instrument has to be worked with intermittently using little pressure and subject to continuous cooling using a sterile physiological saline solution.

ATTENTION: Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Hager & Meisinger GmbH
Hansemanstr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

Meisinger USA, L.L.C.
10150 E. Easter Avenue
Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400
Fax: +1 (303) 268-5407
Email: info@meisingerusa.com
Internet: www.meisingerusa.com

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