## Lift-Control Systeme

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## Lift-Control Systems

### Systems for Internal and External Sinus Lift

Art.-No.: BTL00, CSL00, BLIPL

The Bone Management®-Systems External-Lift-Control, Internal-Lift-Control und Lift-Control Plus offer instruments for the external sinus lift according to Boyne (direct sinus lift) as well as for the internal sinus lift according to Summers (indirect sinus lift).

### Indications for Use

The Bone Management®-Systems are indicated for performing the external and internal sinus lift in cases of insufficient bone height in the maxillary posterior region. In the case of the internal sinus lift, augmentations of small vertical defects of 3 mm to 4 mm are possible with a residual bone height of ≥ 5 mm.

### Utilization

The specific system compositions as well as the recommended and maximum permitted speeds for the application are summarised in the article overview under „Contents“. For proper, safe use, it is essential to observe these.

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In the first step, the gingiva is folded down with a suitable incision in the area of the maxillary sinus so that the maxillary bone is exposed for further treatment. Using the spherical diamond and carbide instruments included in the set, a lateral window to the maxillary sinus is carefully prepared by thinning the bone in this area down to the membrane. The fine mucous membrane of the maxillary sinus (Schneider's membrane) should not be perforated during the process.

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The Schneiderian membrane is dissected off using the appropriate manual instruments. Depending on the anatomical conditions, the manual instruments are available in different sizes, shapes and angles to carefully loosen the membrane without injury and thus lift it cranially.

The previously opened bone window of the maxillary sinus wall can now be carefully moved upwards into the maxillary sinus with the mucous membrane already detached.

#### Filling the cavity with bone substitute material

The cavity created in the maxillary sinus is now loosely filled with bone replacement material according to the manufacturer's instructions. The augmentation should hold the previously lifted bone window.

The TL8 condenser can now be used to carefully compact the inserted augmentation. The still open lateral window is covered with the selected membrane according to the manufacturer's instructions. If necessary, the TL20 manual instrument can be used to position the membrane. The mucoperiosteal flap is then sutured in a saliva-proof and tension-free manner.

#### Implant insertion and further treatment

Depending on the patient's situation, implants can be placed in the same session (one-stage) or after the augmentation has healed (two-stage). The prerequisite for the single-stage procedure is the achievement of primary stability of the implant. This usually requires at least 3 mm of bone height. In the two-stage procedure, the implantation takes place only after complete ossification of the augmentation. After appropriate healing, the implant(s) is/are uncovered and the prosthetic treatment is placed. The corresponding manufacturer's instructions must be observed.

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Through the use of the pilot drill, the subsequently used twist drill can be prevented from slipping, especially with a pointed alveolar ridge.

#### Preparation of the implant bed for indirect access to the maxillary sinus

The use of the pilot drilling can prevent the subsequently used twist drill from slipping, especially with a pointed alveolar ridge.

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#### Compaction of the jawbone and raising of the sinus floor

In the last step of the preparation, the bone is first compacted laterally and apically with the help of the corresponding elevator. The compacted bone is then driven cranially into the sinus cavity. At the same time, the Schneiderian membrane is raised, whereby care must be taken not to damage the membrane. Driving in the elevators should be done gently and by means of controlled hammer blows. This should be done by means of many short and gentle rather than less forceful hammer blows. If there is a strong cortical bone in the crestal area, it is recommended to widen it with the help of the pilot bore and thus facilitate the preparation with the elevators.

Depending on the individual patient situation, the choice of instruments used may vary. The following applies:

Soft bone / condensation desired:  
maximum drill diameter < maximum elevator diameter

Hard bone / condensation not desired:  
maximum drill diameter > maximum elevator diameter

The diameter refers to the drill or elevator used in the previous preparation step.

In order to achieve an even higher compaction of the bone, if required, it is possible to prepare a smaller diameter hole and only achieve the desired diameter by using several elevators with increasing diameters. The preparation with the elevators is first carried out up to about 1 mm below the sinus floor. Only with the last elevator inserted is the sinus floor, including the bone compacted apically, driven cranially into the sinus cavity, with a simultaneous elevation of Schneider's membrane.

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#### Preparation of the Schneiderian membrane

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EN

**Implant insertion****Possibility 1: Single-stage procedure**

In the single-stage procedure, the newly created space is filled with the selected bone substitute material according to the manufacturer's instructions. Subsequently, the widening and fine preparation of the implant bed can be carried out with the system's own instruments, as well as the insertion of the implant according to the manufacturer's instructions.

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EN

**Possibility 2: Two-stage procedure**

In the two-stage procedure, the first step is to fill the newly created space with the selected bone substitute material according to the manufacturer's instructions. The implantation is carried out according to the manufacturer's instructions in a second procedure as soon as the new bone formation is complete.

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**Hager & Meisinger GmbH**