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Indications for Use

The Bone Management® Transfer-Control Systems are designed for the oral harvesting of bone cylinders and bone rings for the augmentation of both horizontal and vertical bone deficits in the upper and lower jaw.


Utilization

The specific system compositions as well as the recommended and maximum permitted speeds for the application are summarised in the article overview under „Contents“. For proper, safe use, it is essential to observe these.

Harvesting and transplantation of bone cylinders (CTRBA, CTRB2, CTR00, CTRPL)


Horizontal augmentation
In advanced horizontal bone resorption - when bone spreading and bone splitting are no longer possible - a **horizontal augmentation** is required to compensate.


Extraction of the bone cylinder
Depending on the size of the defect and the required implant diameter, a diameter-congruent autologous bone cylinder is first harvested with a suitable trephine, such as from the retromolar region of the lower jaw. If necessary, this can be detached from the donor bone with one of the diamond-coated saws. If necessary, the donor region can be filled with collagen or bone replacement material after bone harvesting according to the manufacturer's instructions. A conventional suture closure ends the procedure.

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

Preparation of the bone site
The bone site on the horizontally resorbed alveolar process is now pre-milled in the tooth axis with the ablative bur matching the previously used trephine. The bur is designed in such a way that the recipient site is well refreshed, thus creating favourable conditions for healing.

Transplantation of the bone cylinder
The harvested bone cylinder from the donor region is now inserted into the recipient site to compensate the horizontal bone deficit. Care must be taken to ensure that the bone cylinder is inserted quickly. Intermediate storage of the transplant takes place in blood or isotonic saline solution.

Fixation of the bone cylinder
Depending on the case situation, the bone cylinder can be fixed with a fixation screw, such as from the Screw System TX (Art.-Nr. BTX00, BTXPR). The manufacturer's instructions for the selected screw system must be followed. If a secure fixation of the transplant with one screw is not possible, fixation with at least 2 screws is recommended to prevent any rotation. If necessary, transitions between the bone site and the transplant can be filled with bone replacement material according to the manufacturer's instructions. The wound closure should be complete and tension-free.


Vertical augmentation
If vertical bone resorption is advanced, a **vertical augmentation** is required to compensate. The form-congruent recipient site is prepared on the crestal alveolar process with the corresponding wheel cutter. Here, too, care must be taken to use instruments with congruent diameters. The harvesting as well as the fixation of the bone cylinder is carried out according to the same principle as for the horizontal augmentation.


Harvesting and transplantation of bone rings (CTRIO, CTRPL)

Pilot drilling of the donor region
First, the donor region is marked with the help of the pilot drill. The surface of the future bone ring is milled flat and a ring groove is created. The pilot drilling achieved with the help of the mandrel is later used to precisely centre the osteosynthesis screw in the transplant.

Extraction of the bone ring
A bone ring of congruent diameter is then milled into the donor region with the aid of the trephine. The depth of the ring depends on the vertical build-ups required.

With the aid of a manual instrument (e.g. from the Surgical Kit BCS00), the bone ring is now carefully removed from the donor region. If necessary, the donor region can be filled with collagen or bone replacement material after bone harvesting according to the manufacturer's instructions. A conventional suture closure ends the procedure. The removed bone should be used as soon as possible. Intermediate storage of the transplant takes place in blood or isotonic saline solution.

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


**Preparing the recipient site**

The bone site on the vertically resorbed alveolar process is now prepared with the ablative bur matching the previously used trephine in such a way that a flat and refreshed contact surface is created, which creates favourable conditions for rapid healing. Thanks to the harmonised instruments, an optimal fit of the transplant is guaranteed.

**Transplantation and fixation of the bone ring**

After inserting the bone ring into the recipient site, it can be fixed with an osteosynthesis screw, such as from the Screw System TX (Art.-No. BTX00, BTXPR).
The manufacturer's instructions for the selected screw system must be followed. If necessary, transitions between the bone site and the transplant can be filled with bone replacement material according to the manufacturer's instructions. The wound closure should be complete and tension-free.

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