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Application and safety instructions for the MEISINGER

**EN**

## Universal torque ratchet

10-70 Ncm

Art.-No.: CARAS


### Indications for Use

Torque ratchet for inserting and removing dental screws with a defined torque. The torque function can also be "blocked"; the blocked position enables greater torque to be transferred when placing implants, and allows connections to be loosened. The torque ratchet may only be used by trained dental specialists.

### Utilization


#### Prosthesis adjustment – function of torque

Application by using the bending rod. The torque will be read by means of the bending rod on the scale. The requested torque will be achieved when the middle of the bending rod will be covered with the appropriate scale graduation mark. ATTENTION: Reading should always be done directly from the top. When the desired torque has been achieved, please relieve the bending rod again. Then, the bending rod will spring back into the starting position.

#### Surgical adjustment – blocked function

Use the ratchet without the bending rod. ATTENTION: The ratchet should not be charged over 100 Ncm.

**Switching over of the direction of rotation:**

- Pull the turning grip
- Turn half way the turning grip
- Release the turning grip

#### Exchanging of the tools

Pull the turning grip, then the tools can be inserted or removed. Then release the turning grip. Tools can be used from both sides. If necessary, switch over the direction of rotation.

### Preparation

#### Treatment instructions / warnings

To avoid damage, do not use metal brushes or cleaning sponges. Only use cleaning and disinfectant solutions with a pH value of between 4.5 and 10. Follow the manufacturer's instructions (e.g. intended purpose, dosage, exposure period and replacement of the solution).


The ratchet is not sterile when delivered and must be cleaned and sterilised before it is used. When using several torque ratchets, do not interchange the individual parts. Each individual part belongs to the respective instrument. Damaged products must go through the reprocessing process before being returned for repair.

#### Restriction regarding reprocessing

The end of the product's service life is normally determined by wear and damage caused during use and by incorrect handling.

#### Preparations for cleaning

For cleaning the torque ratchet can be disassembled without any tool. Please unscrew the screw completely. Then, the entire pawl can be taken out.

Clean the parts under cold running water using a soft brush to remove all visible soiling. Ensure that all openings and cavities are thoroughly rinsed. Do not allow blood and other soiling to dry on.

#### Cleaning and disinfection:

Manual:  
Ultrasonic cleaning bath: Place the parts into a wire basket and ensure that the parts do not touch, in order to avoid acoustic shadows. Clean for 3 minutes in an ultrasonic cleaning bath (35-40 kHz) at a temperature of 40°-50°C with an enzymatic cleaning solution. Ensure that the parts are completely immersed in the water, without the formation of any bubbles. Rinse with clear, cold water; if possible, use deionised water. Dry the parts with a lint-free cloth and blow dry with compressed air.

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